



IC333 Bed

The IC333 Bed is unmatched for function and durability, taking the hospital look out of your room. This four function bed, features dual lift motors totalling 12,000 newtons lifting power, unlimited under-bed clearance and eight heavy duty castors.



KEY FUNCTIONS:



Head Lift



Knee Break



Height Adjustable
22-66cm



Trendelenburg and
Reverse Trendelenburg

BED CODES & SIZES:



LONG SINGLE

Code: IC333LS
Mattress Surface Size: 2030 x 900
Base Size: 2150 x 960



LONG DOUBLE

Code: IC333D
Mattress Surface Size: 2030 x 1350
Base Size: 2150 x 1400
One hand control



KING SINGLE

Code: IC333KS
Mattress Surface Size: 2030 x 1070
Base Size: 2150 x 1120



QUEEN

Code: IC333Q
Mattress Surface Size: 2030 x 1500
Base Size: 2150 x 1580
One hand control

Note: Always remember that TWO Long Single beds together make a King size bed.