




BodyAdjust Sleep System



MULTiWAY

As specialists in sleep systems, we believe that a bed is an important aspect of our daily care. It is a place where rest and recovery should be comfortable, therapeutic & secure.

BodyAdjust

Sleep System by **MULTiWAY**

The **BodyAdjust Sleep System** is designed to blend into the home environment with ease, yet not compromise on your therapeutic needs. Multiple positioning options will assist you with independent mobility, improve sleep quality, and comfort levels.

Functions

- Electric Hi-Lo movement
- Electric back raise
- Electric knee break and leg raise
- Trendelenburg & reverse Trendelenburg positioning
- "Zero Gravity" functionality
- Supine positioning

Features

- Back-lit, wipeable corded handset
- Foot and side mattress retainers (incl)
- Companion capability

Sizes available

- Long Single: 2030mm (L) x 915mm (W)
- King Single: 2030mm (L) x 1070mm (W)
- Double: 2030mm (L) x 1370mm (W)
- Queen: 2030mm (L) x 1525mm (W)
- Split King: 2030mm (L) x 1830mm (W)
- Split Super King: 2030mm (L) x 2140mm (W)

Specifications

- Colour: Midnight Black
- Single & King Single models: 140kg SWL
- Double & Queen models: 140kg SWL
- Hi-Lo lift range: 690mm
- Floor to bottom of base: 155 mm (in lowest position)
- Floor to top of base: 240mm (in lowest position)
- Trendelenburg tilt angles: 17 degrees
- Maximum back raise angle: 72 degrees
- Maximum leg raise angle: 46 degrees
- Hoist clearance: 900mm
- Castors: 55mm x 8 (4 x lockable)
- Warranty: 5 x years frame, 5 x years electrical

Accessories (optional)

- Headboard
- Side grip handles
- ¾ length bed rails

Mattress Options

Active

A multi-layered mattress giving you all the benefits of contouring foam, the breathability of cotton and bamboo fibers, and the comfort of pocket springs. Ideal for those requiring minimal partner disturbance through the night.

Comfort

This supple 2-layer combination foam mattress with latex upper supports and conforms to the contours of your body closely. The natural latex allows air to circulate, and the immersion/support of the foam/latex together aids in pressure redistribution.

Mattress Overlay

A beautiful, high-quality 50mm memory foam overlay to immerse into - providing pressure point relief, or longevity to your existing mattress.