



The Entrix Range



Integrated pressure monitoring mattress overlay system for all your essential needs.

Simple but powerful, the Carilex Entrix is a cost effective therapy support surface that offers a combination of dynamic and static modes to meet individual patient's need. Working with the state-of-the-art IPS technology, Entrix's intelligent pressure monitoring system keeps an eye on the mattress 24 hours a day, ensuring the patient receives pressure relief or pressure redistribution at all times.

Adjustable weight settings for different body types, low pressure alert in case of abnormal pressure drop or power failure and Two Directional Stretch Textile for durability and comfort.

	ENTRIX 5" OVERLAY	ENTRIX 7" MATTRESS REPLACEMENT
Risk	Low to Medium	Low to Medium
SWL	150kg	150kg
CPR	✓	✓
Cable Management	✓	✓
Cells	16 x 5" cells	16 x 5" cells with 2" built in foam
Two Way Direction Stretch Textile	✓	✓
Low Pressure and Power Failure Alarms	✓	✓
Cycle Times	10 minutes	10 minutes
Mattress Sizes	Narrow – 82cm x 200cm Single – 90cm x 200cm King Single – 107cm x 200cm	Single – 90cm x 200cm King Single – 107cm x 200cm Euro King Single – 100cm x 200cm
Cost	\$\$	\$\$\$