



## Breezy EC 2000

A quality wheelchair with options offering the flexibility to accommodate a variety of users.

### Standard Features

- Dual axle allows the seat to be easily lowered to hemi height
- Chrome-plated, carbon steel frame is built to last, great for rental fleets
- Composite mag wheels are molded with contours to enhance strength and durability
- Full size chart pocket
- Industry-compatible front riggings
- 250 lb. weight capacity (300 lb. on 20" seat width)

### Options for Versatility

- Removable desk-length or full-length or fixed full-length armrests
- Swing-away footrests, elevating legrests or fixed footrests

### Models

216FAFPFS	16" width, fixed padded full arms, fixed footrests
216FADPS	16" width, fixed padded full arms, swing-away footrests
216RADPS	16" width, removable padded desk arms, swing-awa footrests
216RADPE	16" width, removable padded desk arms, elevating legrests
216RADUS	16" width, removable desk arms, swing-away footrests
216RAFPS	16" width, removable padded full arms, swing-away footrests
216RAFPE	16" width, removable padded full arms, elevating legrests
218FAFPFS	18" width, fixed padded full arms, fixed footrests
218FADPS	18" width, fixed padded full arms, swing-away footrests
218RADPS	18" width, removable padded desk arms, swing-away footrests
218RADPE	18" width, removable padded desk arms, elevating legrests
218RADUS	18" width, removable desk arms, swing-away footrests
218RAFPS	18" width, removable padded full arms, swing-away footrests
218RAFPE	18" width, removable padded full arms, elevating legrests
220RADUS	20" width, desk arms, swing-away footrests
220RADPS	20" width, padded desk arms, swing-away footrests
220RADPE	20" width, padded desk arms, elevating legrests
220RAFPS	20" width, padded full arms, swing-away footrests
220RAFPE	20" width, padded full arms, elevating legrests