



## Karma Ergo Lite Deluxe



At 10.3kg the Karma Ergo Lite Deluxe is the lightest self propelled wheelchair in Karma's range. It uses Karma's S Ergo patented seating system to help prevent the user slipping off the seat and adjusts and stabilises seating pose, with height adjustable, swing-away footrests. Folding, lightweight frame and fixed ergonomic arm rest pad as standard. Available in Transit also.



## Karma S-Ergo Self-Propelled



### Features

- Lightweight aluminium frame
- Swing-away removable footrest
- Swing-away armrest with ergonomic pads
- Quick release solid tyre wheels
- Solid 7" castors
- Push-to-lock with attendant braking system
- Fixed seat height
- Mesh with aegis microbe shield padding upholstery

## Karma S-Ergo Transit



### Features

- Folding Lightweight aluminum frame
- Swing-away removable footrest
- Swing-away with ergonomic pads
- Fixed solid tyre wheels 14" x 1.75"
- Solid 7" castors
- Push-to-lock with attendant braking system
- Fixed seat height
- Mesh with aegis microbe shield padding upholstery