FOLD & GO

COMPACT FOLDING POWERCHAIR P113





www.meritsaustralia.com

FEATURES

- Lightweight, foldable power wheel chair which folds easily in seconds for transportation or storage
- Folds in one easy step
- Compact size means great manoeuvrability in tight spaces
- Removable battery box with battery charge indicator & extra charge port for charging off board
- · Adjustable height footplate
- · Adjustable height backrest
- Intelligent, regenerative, electromagnetic braking
- · Vibrant smart colours available, red, blue or grey

SPECIFICATION

Max. Range	11km / 7mi. 15km / 9.6mi. (optional)	
Max. Speed	7.5kph / 4.7kph	
Battery	1pc Lithium-ion battery 24V/11.4AH / 17.1AH (optional)	
Motor	DC 24V, 70W (Rating)/210W (Max)	
Gradient	3°	
Controller	PG VR2 50A	
Charger	2A	
Caster	7" PU tire	
Drive Wheel	8"(200x50) PU tire	
Brake	electromagnetic brakes	

DIMENSIONS (cm/in)



WEIGHT (kg/lb)

ttl w/batt.	ttl w/o batt.	capacity
23/50.5 23 9/52 5 (optional)	20.9/46	114/250

COMPACT FOLDING POWER WHEELCHAIR







www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.





