

ProCair Plus



homecare equipment

Alternating pressure solutions for the prevention and treatment of pressure injuries.



- ❑ Fully automated operational features
- ❑ Advanced safety features including 5 hour back-up battery
- ❑ Advanced error and service indicators
- ❑ Optimal patient therapy and comfort



FAILSAFE

Inbuilt backup battery to provide normal operation for up to 5 hours



AUTOFWLER

Automatically boost pressure to protect the sacral area during seated position




AUTOCARE

Automatic weight sensing for ease of use and optimal patient comfort



HEELCAIR

7 disconnectable narrower heel cells provide full offloading for heels



SAFETY CARE

Cell in cell construction providing safety from bottoming out

PURE TPU

Design for softer, quieter and more comfortable experience for the patient

CAIRALERT

Bed exit for falls risk management

QUICK RELEASE CPR

Easy opening and accessible CPR unit for rapid deflation under 10 seconds

SPECIFICATIONS

CODE	APMPC-R02S	APMPC-R02SK
SIZE	Single	King Single
CAPACITY (SAFE WORKING LOAD)	40-220 kg	40-220 kg
SYSTEM ALTERNATION	1 in 2	1 in 2
WARRANTY	3 year	3 year
PRESSURE INJURY RISK	High to very high risk	High to very high risk
ARTG	289458	289458
MATTRESS DIMENSIONS (L x W x H)	2000 x 880 x 200 mm	2000 x 1050 x 200 mm
MATTRESS WEIGHT	13 kg	13.5 kg