



HC3 V2

Three Layer Mattress

Features & Benefits

- Shear reducing 'MEMAIR VISCO' castellated top layer
- Firm perimeter side walls to assist with patient transfer
- 'Purify' antimicrobial foam treatment
- Fire retardant foam and cover
- BS6807 & BS7175 - Crib 5 - PASS
- 50mm 60/50 high density 'MEMAIR VISCO' castellated top layer
- Medium density middle layer
- High density base layer
- 1.8 support factor
- Heel slope with visco elastic foam
- Omni-directional stretch top cover with PVC base cover
- 3 sided zip with waterfall covered flap

Specifications

CODE	TYPE	WIDTH	LENGTH	HEIGHT	MAX USER WEIGHT
BDM00120	Single	860mm	2000mm	150mm	200kg
BDM00121	King Single	1050mm	2000mm	150mm	200kg