

SALSA 200

Compact Active Standing Aid



Allegro
CONCEPTS PTY LTD

Innovation in Patient Handling

SALSA 200

SALSA 200: Designed as an active standing aid allowing a single carer to raise a patient into a comfortable secure standing position.

The Salsa 200 is ideal for transferring from bed to chair, or chair to commode. Allegro Concepts standing slings are manufactured to suit the Salsa, allowing superior access for dressing and personal hygiene.

With a small footprint and low overall weight, the Salsa 200 is highly maneuverable. A high capacity safe working load of 200kg means the Salsa is equally at home in an institutional or home environment.

Salsa 200 Specifications

| | |
|----------------------------|-------------------|
| Safe Working Load | 200kg |
| Base Length | 1010mm |
| External Base Width Closed | 665mm |
| External Base Width Open | 1030mm |
| Internal Base Width Closed | 565mm |
| Internal Base Width Open | 930mm |
| Base Width Adjustment | Electric |
| Under Bed Clearance | 130mm |
| Minimum Height Lift Arm | 880mm |
| Maximum Height Lift Arm | 1685mm |
| Under Base Height | 55mm |
| Lifter Weight | 51kg |
| Minimum Storage Height | 1190mm |
| Footplate | Removable |
| Kneepad | Height Adjustable |



Multiple Grip Positions




Adjustable Kneepad



Removable Footplate



Tempo Handle Kit - THK01



Custom Sitting Sling



Knob Kit for Custom Sitting Sling - KK01

The Allegro Concepts range of patient lifting equipment including lifters, slings and accessories has been designed and tested to comply with: AS ISO 10535-2011 Hoists for the transfer of disabled persons – Requirements and test methods.

Allegro Concepts Pty Ltd

Unit 9, 110- 120 Silverwater Road | Silverwater | NSW | 2128

☎ 02 6895 2888 | ✉ sales@allegroconcepts.com.au

🌐 www.allegroconcepts.com.au