



BodyAdjust Sleep System by MULTiWAY

About MULTiWAY

MULTiWAY is an **Australian-owned** company that prides itself on providing quality imported products at competitive prices.


We work closely with healthcare and equipment providers, prescribers and end-users, to identify the latest technologies and changing needs of the market and our clients.

As specialists in sleep systems, we believe that a bed is an important aspect of our daily care. It is a place where rest and recovery should be comfortable, therapeutic and secure.

The **BodyAdjust Sleep System** is designed to blend into the home environment with ease, yet not compromise on your therapeutic needs. Multiple positioning options will assist you with independent mobility, improve sleep quality, and comfort levels.

Voice ACTIVATION

Sizes Available



Mattress Options

Active


A multi-layered mattress giving you all the benefits of contouring foam, the breathability of cotton and bamboo fibres, and the comfort of pocket-springs. Ideal for those requiring minimal partner disturbance through the night.

Comfort

This supple 2-layer combination foam mattress supports and conforms to the contours of your body closely. The natural latex upper-layer allows air to circulate and improves immersion to assist with pressure redistribution.

Mattress Overlay


A beautiful, high quality 50mm memory foam overlay to immerse into - providing additional comfort to your existing mattress.



with Smartphone App Control

BodyAdjust

Sleep System by MULTiWAY



SPECIFICATIONS

Colour	Midnight Black
Safe Working Load	Single & King Single models: 195kg Double & Queen models: 195kg
Hi-Lo Lift range	690mm
Floor to bottom of base	155mm (in lowest position)
Floor to top of base	240mm (in lowest position)
Trendelenburg tilt angles	17 degrees
Maximum back raise angle	72 degrees
Maximum leg raise angle	46 degrees
Hoist clearance	Yes
Castors	55mm x 8 (4 lockable)
Warranty	5 year - Frame 2 year - Electrical



MULTiWAY

CONTACT US

P. 1300 168 584 (1300 1 MULTI) M. 0411 730 124

E. sales@multiway.com.au W. multiway.com.au