

# Ultimate Comfort

## CURA CONTOUR POSITIONING CHAIR

Experience unparalleled comfort and support with the Cura Contour, the world's most adjustable positioning chair. Designed to address even the most complex postural needs, this chair offers exceptional flexibility for complex postural support.

Revolutionary  
25-point Adjustable  
Backrest System



homecare equipment

Homecare Equipment

P: 08 8338 7988

E: [info@homecareequipment.com.au](mailto:info@homecareequipment.com.au)

[homecareequipment.com.au](http://homecareequipment.com.au)

# Cura Contour Chair



## PRODUCT OVERVIEW

The Cura Contour chair sets a new standard in adjustability with its Oyster backrest system, offering 25 unique adjustment points to easily accommodate severe postural deformities. Available in two sizes for growth adaptability, it's a smart investment. This range comes in attractive standard colors, with extensive fabric and color options for special orders.

## CHAIR OPTIONS

- ✓ Fully adjustable system
- ✓ Electric operation
- ✓ Padded footplate
- ✓ Attractive colour options
- ✓ Optional cushion upgrades


Model	Kinder	Adult
Parameter	Inch/cm	Inch/cm
Seat Width Range	12/30.5 to 18/45.7	16/40.6 to 22/55
Seat Depth Range	12/30.5 to 20/50.8	16/40.6 to 22/56
Backrest Height	26/66	30/76.2
Seat to Footplate	8/20.3 to 25/63.6	8/20.3 to 25/63.6
Seat to Armrest Height	4/10 to 6/15	7/18 to 10/25
Tilt in Space Range	0 to 30 Degrees	0 to 30 Degrees
Seat to Back Angle Range	-5 to 40 Degrees	-5 to 40 Degrees
Floor to Seat Height	19/48.3 to 22/56	16/40.6 to 22/56
Maximum User Weight	160kg	160kg

## SUITABLE FOR:

MS, MND, Paraplegia, Scoliosis, CP, Lordosis, Kyphosis, Elderly, Neurological Debilitating Conditions.

## The World's Most Adjustable Positioning Chair

Experience unmatched versatility with the Cura Contour Chair, equipped with a 25-point adjustable Oyster backrest, designed to meet extreme postural challenges and adapt to growth and evolving needs.



## Homecare Equipment

P: 08 8338 7988

E: [info@homecareequipment.com.au](mailto:info@homecareequipment.com.au)

[homecareequipment.com.au](http://homecareequipment.com.au)

