

Molift Duo

A flexible and free standing gantry system for temporary hoist situations

The easy solution for most room sizes

Molift Duo offers a system that is easy to install without being fixed to the ceiling or walls. It is ideally suited to temporary hoist applications in hospitals, nursing homes or private homes. In combination with Molift Air or Molift Nomad hoist and a sling, Molift Duo solves most overhead lifting situations, whatever the environment. With a lifting capacity of 300 kg, it is ideal for lifting and transferring a person from or to; bed, floor, chair, wheelchair, toilet etc.

Lightweight construction


Molift Duo has a lightweight aluminium construction and a functional design which makes the gantry system easy to install, dismantle and move elsewhere.

This 2-post gantry system has a rail length between 2000 mm and 3500 mm and has 4 step adjustable heights up to 2550 mm, making this free standing hoist unit suitable for all bed and room sizes.

Features

- Easy to install - due to the lightweight system, it is easy to install the gantry and move it to a new location.
- Supreme stability - the gantry system has supreme stability and a functional design.
- Adjustable height - height positions can be set at four fixed positions from 225 cm up to 255 cm. Adjust the height with the quick release pin.
- Charging of lift motor - Molift hoist motors can easily be charged through the hand control or optionally through in-rail charging.



Description	Item no.	Dimensions	Material
Molift Duo 2000	1600013	Rail lengths 2000/2500/3000/3500 mm	Anodized aluminium, Plastic  For more information visit etac.com
Molift Duo 2500	1600015	Height 2250 - 2550 mm	
Molift Duo 3000	1600017	Depth 1200 mm	
Molift Duo 3500	1600019	Maximum internal space 3500 x 2400 mm	
		Weight Total weight (complete unit) 38 - 44,3 kg	

